

Physioage Medical Group
 30 CENTRAL PARK S, SUITE 8D
 NEW YORK NY 10019

NAME: Sample Report	AGE:
External Client: Physioage Medical Group	SEX: Unknown
External Client MRN: ...	DOB:
	UCLA MRN: N/A

Authorizing Provider:

Copy To:

See Values: T-Helper/Inducer (CD4+) (L), CD8+/CD95- in CD3 (CD3+CD8+CD95-) (L), CD8+/CD95- in CD8 (CD8+CD95-) (L)

ID: 24IM-004100203
 Collected: 1/3/2024 1030
 Received: 1/4/2024 1333
 Verified On: 1/9/2024 1449

Type/Src: Blood
 Authorized by: Raffaele, Joseph, MD
 Resulting Lab: UCLAImmuno

Senescent T CD8 (Final result)

Component	Value	Ref. Range
T-Lymphocytes (CD3+)	67	49 - 84 %
T-Lymphocytes (CD3+)	775	603 - 2,990 /cmm
T-Helper/Inducer (CD4+)	35	28 - 63 %
T-Helper/Inducer (CD4+)	406	(L) 441 - 2,156 /cmm
T-Suppressor/Cytotoxic (CD8+)	29	10 - 40 %
T-Suppressor/Cytotoxic (CD8+)	330	125 - 1,312 /cmm
Ratio (CD4:CD8)	1.23	0.96 - 3.93
B-Lymphocytes (CD19+)	20	7 - 27 %
B-Lymphocytes (CD19+)	225	107 - 698 /cmm
NK (CD56+/CD16+)	12	4 - 25 %
NK (CD56+/CD16+)	142	95 - 640 /cmm
CD8+/CD28- in CD3 (CD3+CD8+CD28-)	19	1 - 28 %
CD8+/CD28- in CD3 (CD3+CD8+CD28-)	147	11 - 359 /cmm
CD8+/CD95- in CD3 (CD3+CD8+CD95-)	4	3 - 27 %
CD8+/CD95- in CD3 (CD3+CD8+CD95-)	31	(L) 33 - 357 /cmm
CD8+/CD28- in CD8 (CD8+CD28-)	48	4 - 51 %
CD8+/CD28- in CD8 (CD8+CD28-)	158	17 - 364 /cmm
CD8+/CD95- in CD8 (CD8+CD95-)	10	(L) 11 - 57 %
CD8+/CD95- in CD8 (CD8+CD95-)	33	32 - 347 /cmm

Test Comment
 Testing performed by flow cytometry.